

RACCOON

Steam Sterilizer

(Terminal Steam Sterilizer)

MHS-70EV

■ *Sterilizing milk preparation utensils*

■ *Terminal sterilizing and cooling baby milk also*



<< SPECIFICATIONS >>

Dimension	W 700mm x D 800mm x H 1000mm
Material	Body : Stainless steel SUS-304
Capacity	4 containers or 4 baskets (24 standing feeding bottles each)
Sterilization purpose	1. Terminal sterilization mode : Terminal sterilization for baby milk 2. Utensil sterilization mode : Milk preparation utensils including nipples and caps
Control	Sterilization by high temperature water
Temperature control	By heat transfer sensor
Sterilization temperature	Settable by 1degree (max 95 degrees Celsius)
Sterilization time	Settable by second
Safety system	Ejector steam termination system Boil-dry protection
Lid	Stainless steel, lock mechanism while open
Plumbing	Cold water supply : 20A Hot water supply : 20A Drainage : 40A (High temperature drainage)
Power source	AC200V three-phase 4Kw (50Hz/60Hz)

RACCOON

MITA RIKA KOGYO CO., LTD.

2-8-2, Oyodo-naka, Kita-ku, Osaka, 531-0076, Japan

Tel : 81-6-6458-0971

E-mail : info@racoon.co.jp

Fax : 81-6-6458-2393

URL : www.racoon.co.jp

