



Model DW-FL90 / DW-FL135 / DW-FL270

Features

- Vertical cabinet (upright)
- Operating temperature ranging from -10°C to -40°C can be set freely
- Microprocessor temperature controller
- LCD display
- Door with lock
- Audible / visual alarm : High or low temperature alarm, system failure alarm
- Danfoss compressor and EBM fan motor
- Password protected configuration page (model DW-FL262 and DW-FL362 only)
- 4 casters provided for easy mobility
- Refrigerant : R507, CFC free
- Certificates : CE mark for DW-FL 90/270, ISO9001, ISO14001, ISO13485
- OPTIONAL : 7- day chart recorder



Model DW-FL262 / DW-FL362

Specifications

Model	DW-FL90	DW-FL135	DW-FL270	DW-FL262	DW-FL362
CE MARK	Yes	-	Yes	-	-
Capacity	90 litres	135 litres	270 litres	262 litres	362 litres
Drawers/shelves (fixed position)	4 drawers	5 drawers	7 drawers	4 shelves	5 shelves
Exterior construction	Powder coated mild – steel exterior				
Cabinet construction	ABS interior			Stainless steel interior except door	
Insulation	90 mm polyurethane foam			110 mm polyurethane foam	
External size mm	623W x 544D x 842H	623W x 544D x 1123H	623W x 544D x 1680H	930W x 731D x 1227H	930W x 731D x 1627H
Internal size mm	438W x 388D x 539H	438W x 388D x 794H	438W x 388D x 1580H	580W x 426D x 1036H	658W x 426D x 1290H
Input power	280 watt	280 watt	370 watt	300 watt	400 watt
Electrical	220 – 240 VAC, 50 Hz, single phase				
Weight	44/48 kg	54/60 kg	65/70 kg	86/92 kg	102/110 kg